

HANAROSTENT®

TwynLay® (NCN)



Indication :

- Palliative treatment of pyloric or duodenal stricture caused by malignant tumors

Features :

- Flare heads on the both ends for Anti-Migration
- Uncovered on the both ends for Anti-Migration
- Double layered body structure and high conformity

Fluoroscopic Application


	Stent (mm)			Delivery device	
	Diameter	Usable length	Total length*	Length(mm)	Diameter(mm/fr)
TLD-20-060-180	23-20-23	30	60	1800	3.4/10.2
TLD-20-170-180		140	170	1800	3.4/10.2
TLD-22-060-180	26-22-26	30	60	1800	3.4/10.2
TLD-22-170-180		140	170	1800	3.4/10.2
TLD-24-060-180	30-24-30	30	60	1800	3.4/10.2
TLD-24-150-180		120	150	1800	3.4/10.2

*Increment of 10mm

Fluoroscopic Application


	Stent (mm)			Delivery device	
	Diameter	Usable length	Total length*	Length(mm)	Diameter(mm/fr)
TLD-20-060-120	23-20-23	30	60	1200	3.4/10.2
TLD-20-170-120		140	170	1200	3.4/10.2
TLD-22-060-120	26-22-26	30	60	1200	3.4/10.2
TLD-22-170-120		140	170	1200	3.4/10.2
TLD-24-060-120	30-24-30	30	60	1200	3.4/10.2
TLD-24-150-120		120	150	1200	3.4/10.2

*Increment of 10mm

Distal  Proximal	(NNN) Non covered	<u>XXX-00-000-000</u>	<u>XXX 00 000-X000</u>
	(CCC) Fully covered	① ② ③ ④	① ② ③ ④
	(NCN) Partially covered	① Product Name	③ Stent Length
		② Stent Diameter	④ Delivery Device Length

HANAROSTENT®

TwynLay® (NC)



Indication :

- Palliative treatment of pyloric or duodenal stricture caused by malignant tumors


Features :

- Flare heads on the both ends for Anti-Migration
- Uncovered on the both ends for Anti-Migration
- Double layered body structure and high conformity

Endoscopic Application

	Stent (mm)			Delivery device	
	Diameter	Usable length	Total length*	Length(mm)	Diameter(mm/fr)
TLD-20-060-230	23-20-23	30	60	2300	3.4/10.2
TLD-20-170-230		120	150	2300	3.4/10.2
TLD-22-060-230	26-22-26	30	60	2300	3.4/10.2
TLD-22-170-230		120	150	2300	3.4/10.2
TLD-24-060-230	30-24-30	30	60	2300	3.4/10.2
TLD-24-150-230		120	150	2300	3.4/10.2

*Increment of 10mm

Distal		Proximal	(NNN) Non covered	<u>XXX-00-000-000</u>	<u>XXX 00 000-X000</u>
			(CCC) Fully covered	① ② ③ ④	① ② ③ ④
			(NCN) Partially covered	① Product Name	③ Stent Length
				② Stent Diameter	④ Delivery Device Length