




SU HANDLE

SY HANDLE

RD HANDLE

RY HANDLE



**STERILE EO**

**Coding**

OSK   5   20   12   SU  
 1        2        3        4        5

- 1- **OSK:** Snare Model
- 2- **5:** Catheter Diameter
- 3- **20:** Snare Opening
- 4- **12:** Catheter Length
- 5- **SU:** Handle Feature

Snare Model	OSK HSK ASK OSW OSB OSWB	Oval Snare Hexagonal Snare Asymmetric Snare Oval Wire Snare Oval Adjustable Snare Oval Wire Adjustable Snare
Catheter Diameter	5 7	5 F (1,8 mm) 7 F (2,3 mm)
Snare Opening	20 25 30 35 40 45 20/30	20 mm 25 mm 30 mm 35 mm 40 mm 45 mm Adjustable Between 20-30 mm
Catheter Length	12 18 21 23 25 31	120 cm 180 cm 215 cm 230 cm 250 cm 310 cm
Handle Feature	SU SY RD RY	Fixed End Fixed End Flushing Channel Rotatable Rotatable Flushing Channel